INSTRUCTIONS FOR USE
PRODUCT SPECIFIC INFORMATION
ONLY ON THIS PAGE

TEGERA® 858

Disposable glove, 0.15 mm nitrile, non powder, Cat. III, purple, extra long, phthalate-free, accelerator-free, latex-free, for touchscreen, waterproof, for precision work

EN 420:2003 + A1:2009  EN ISO 374-1:2016/Type B JKPT

EN ISO 374-5:2016  **VIRUS**

SIZE RANGE (EU) 6,7,8,9,10,11

EU-TYPE EXAMINATION 2777 Satra Technology Europe Ltd Bracetown Business Park, Clonoe, Dublin 15, Dublin, Ireland

ONGOING CONFORMITY CARRIED OUT BY 2777 Satra Technology Europe Ltd Bracetown Business Park, Clonoe, Dublin 15, Dublin, Ireland

TEST ACCORDING TO EN ISO 374-1:2016/ EN ISO 374-4:2019

Tested chemical	Permeation level	Degradation %
J: N-HEPTANE (CAS NUMBER 142-85-5)	2	39,3
K: SODIUM HYDROXIDE 40% (CAS NUMBER 1310-73-2)	6	-1,7
P: HYDROGEN PEROXIDE 30% (CAS NUMBER 772-2-84-1)	3	29,7
T: FORMALDEHYDE 37% (CAS NUMBER 50-00-0)	4	33,4

INSTRUCTIONS FOR USE - CATEGORY III **EN**
 SEE FRONT PAGE FOR PRODUCT SPECIFIC INFORMATION

Carefully read these instructions before using this product.

DECLARATION OF CONFORMITY
 www.ejendals.com/conformity

EXPLANATION OF PICTOGRAMS 0 = Below the minimum performance level for the given individual hazard X= Not submitted to the test or test method not suitable for the glove design or material

Warning! This product is designed to provide protection specified in PPE Regulation (EU) 2016/425 with the detailed levels of performance presented below. However, always remember that no item of PPE can provide full protection and caution must always be taken when exposed to hazardous chemicals or other high risk situations. The performance levels are for products in new condition and do not reflect the actual duration of protection in the workplace due to other factors influencing the performance such as temperature, abrasion, degradation, etc.

EN ISO 374-1:2016 TYPE A, B, C	Protective gloves against dangerous chemicals and microorganisms - Part 1: Terminology and performance requirements for chemical risks. EN ISO 374-1:2016. Definition of break time through the glove palm (Lugpr) (cm ² /min). Type A > level 2 for 6 Chemicals, Type B > level 2 for 3 Chemicals, Type C > level 1 for 1 Chemical.	A: Methanol B: Acetone C: Acetonitril D: Dichloromethan E: Carbon disulfide F: Toluene G: Koltsulfid H: Tetrahydrofuran I: Ethyl acetate	J: n-Heptane K: Sodium hydroxide 40% L: Sulfuric acid 96% M: Nitric acid 65% N: Acetic acid 99% O: Ammoniumhydroxid 25% P: Hydrogen peroxide 30% Q: Fluorsvavelsyra 40% R: Formaldehyd 37%
Permeation level	1 2 3 4 5 6		
Minimum break through times (min)	>10 >30 >60 >120 >240 >480		

ABDEFGH JKLNOPST

Warning: EN ISO 374-1:2016 This information does not reflect the actual duration of protection in the workplace or the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only with respect to the chemical tested. It can be different if used in a mixture. It is recommended to check that the gloves are suitable for the intended use since the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by contact with the chemical, etc., may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider when choosing chemical resistant gloves. Before use, inspect the gloves for any defect or imperfections. For single use only. Degradation is the percentage change in puncture resistance measured after continuous contact with challenge chemical. EN ISO 374-4:2019

EN ISO 374-5:2016 Protective gloves against dangerous chemicals and microorganisms - Part 5 Terminology and performance requirements for microorganism risks. Protection against virus, bacteria and fungi - Pass
Warning: EN ISO 374-5:2016 The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

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Ejendals Ltd, Sweden House, 5 Upper Montagu Street, London, England, W1 2AG
EJENDALS AB